
Examination gloves Nitrile SENSE LIGHT, white

REF code	Product size	PCS/Carton
100018	XS	1600
100019	S	1600
100020	M	1600
100021	L	1600
100022	XL	1440

Resistance to Permeation by Chemicals	Performance Level (EN ISO 374-1:2016 +A1:2018)	Mean Degradation (%) (EN ISO 374-4:2019)
40% Sodium Hydroxide (K)	6	-9.5
37% Formaldehyde (T)	4	16.1

Performance Level	AQL Level	Inspection Level
Level 3	< 0.65	G1
Level 2	< 1.5	G1
Level 1	< 4.0	S4

Permeation Performance Level	Measured Break through time
Level 1	> 10 minutes
Level 2	> 30 minutes
Level 3	> 60 minutes
Level 4	> 120 minutes
Level 5	> 240 minutes
Level 6	> 480 minutes

EN ISO 374-5:2016	EN ISO 374-1:2016 + A1:2018 Type C
VIRUS	
Water tight single use gloves (Personal Protective Equipment) EN ISO 374-5:2016	
- Protection against bacteria and fungus - Pass	
- Protection against viruses - Pass	
These gloves are Personal Protective Equipment of Complex Design Category III, in compliance with Regulation (EU) 2016/425 Type C tested to:	
EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-2:2019, EN ISO 374-4:2019, EN ISO 374-5:2016 and EN 16523-1:2015. The gloves are to act as barrier to protect users against chemicals.	
Medical Device acc. to Regulation (EU) 2017/745.	
Manufactured in accordance with Class I.	
Tested according to EN 455 - 1, 2, 3, & 4.	

en

The penetration resistance has been tested under laboratory conditions, results are related only to the specimen tested. The levels of degradation indicate the change in puncture resistance after exposure to the chemical tested. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical permeation levels are based on tests performed in laboratory using fixed conditions and samples taken from the palm. Result relates only to the chemical tested, the level may vary if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

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