

Gloves



May, 2025

evercare[®]
medical | Part of
Asker
Healthcare
Group

We provide medical supplies to the healthcare sector that make a difference

Evercare Medical develops products that make a positive difference in the European healthcare sector, providing a comprehensive offering within medical supplies to health care providers and patients. We aim to offer a best-in-class portfolio in terms of quality and sustainability, including the new brand Embra® with products that have a lighter environmental footprint, evercare®, and SELEFA®. With over 2,800 products in 33 different product groups, our mission is to support the healthcare sector to improve patient outcomes and reduce total cost of care, while working with our suppliers to ensure sustainable practices throughout the value chain. Our suppliers and manufacturers are carefully selected to ensure they share our view and meet our high requirements in terms of quality, the environment and social responsibility. All products undergo a comprehensive quality control process with extensive quality reviews, including regular inspections of production facilities and randomized audits. Evercare Medical operates two business units, Surgery, Wound & Hospital Care and Gloves & Protection, currently serving 14 countries across Europe. We are headquartered in Stockholm, Sweden and count 35 employees. Evercare Medical is a part of Asker Healthcare Group.

Evercare Medical – we care more.

Embra 
Caring for the future

Embra, derived from the word “embrace” stands for a clear commitment to embrace both environmental and social sustainability. Embra works consistently and continuously to give the world healthcare solutions that focus on lowering environmental impact, while delivering quality and performance for patients and caregivers.

Our mission is to develop responsible healthcare solutions focused on lowering the environmental impact while heightening quality for patients and caregivers. We always consider utilizing innovative processes and raw materials that enable us to reach our goals. We also hold each and every Embra product to strong ethical codes, and to quality reviews throughout its entire value chain.

This is the Embra way. It's our way of caring for patients, caregivers and the future.

 **evercare**
Solutions for premium care

The evercare brand is composed of a premium product range that delivers a comprehensive assortment of examination gloves, surgical gowns and drapes, and wound care, infusion and injection solutions. Offering premium value, more features, and high quality across the product portfolio, evercare is driven by continuous development and innovation to treat patients with the best possible clinical solutions.

 **SELEFA**
Affordable solutions for care

The SELEFA brand provides value for money through its wide range of affordable products, including examination gloves, paper and plastic consumables, non-woven products, bandages and gauze.

Development under the brands Embra®, evercare® and SELEFA® follows a quality management system with established working methods and processes according to the Medical Device Regulation (MDR).

We are powered by a commitment to meeting customer needs, undertaking extensive quality control, and investing in efforts to reduce our impact on the environment, all in the name of improved patient outcomes.



Embra 

 **evercare**

 **SELEFA**

Stay updated




Access our **web page** to find your nearest distributor.



News & Blogs articles with knowledge sharing.




Find the right glove for the right purpose through the **glove guide**.



Watch our **How to use & Donning** guides with videos and instructions.



Find informational videos on **YouTube!**



Follow our **LinkedIn** page and stay updated with our latest news.

Content Gloves




Nitrile Exam Gloves	6
Embra® Nitrile Exam Gloves	7
evercare® Nitrile Exam Gloves	10
SELEFA® Nitrile Exam Gloves	18




Vinyl Exam Gloves.....	24
evercare® Vinyl Exam Gloves	25
SELEFA® Vinyl Exam Gloves	26



Embra® Latex Exam Glove	28
-------------------------------	----



SELEFA® Copolymer Glove	30
-------------------------------	----



SELEFA® TPE Glove	32
-------------------------	----

Nitrile Exam Gloves



Embra® Examination Gloves, Nitrile Proceed Acc Free

REF: 711000, 711001, 711002, 711003, 711004

Nitrile examination gloves and protective gloves for single use, accelerator-free, green.






Key Features

- ▶ Life Cycle Analysis conducted in accordance with the greenhouse gas (GHG) protocols
- ▶ Decreases the total energy usage in the production phase by 55%, based on objective third party tests and LCA reviews, in comparison with Evercare Medical SELEFA® Sense Blue Examination gloves
- ▶ Provides 27% lower CO2 per glove in comparison with SELEFA® Sense Blue Examination gloves, based on a life cycle analysis (LCA)
- ▶ High-performance, lightweight gloves that help reduce waste weight
- ▶ AQL 1.0, ensuring high quality
- ▶ Accelerator-free gloves with low levels of chemical residues, suitable for users with type IV allergy
- ▶ Latex-free and therefore suitable for those with latex / type I allergy
- ▶ Powder-free: reduces risk of allergy and minimises particulate contamination
- ▶ Suitable for contact with all type of foods
- ▶ Light, thin and strong



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
711000	Examination gloves, Nitrile PROCEED, XS, Green		240 mm	200 / 2000
711001	Examination gloves, Nitrile PROCEED, S, Green		240 mm	200 / 2000
711002	Examination gloves, Nitrile PROCEED, M, Green		240 mm	200 / 2000
711003	Examination gloves, Nitrile PROCEED, L, Green		240 mm	200 / 2000
711004	Examination gloves, Nitrile PROCEED, XL, Green		240 mm	180 / 1800

Sales Unit: Transport Carton



Embra® Examination Gloves, Nitrile Proceed X-Long Acc Free

REF: 711005, 711006, 711007, 711008, 711009

Nitrile examination gloves and protective gloves for single use, accelerator-free, green.

Key Features

- ▶ Life Cycle Analysis conducted in accordance with the greenhouse gas (GHG) protocols
- ▶ Decreases the total energy usage in the production phase by 55%, based on objective third party tests and LCA reviews, in comparison between Embra Proceed Examination glove, standard length, with Evercare Medical SELEFA® Sense Blue Examination gloves
- ▶ Provides 27% lower CO2 per glove in comparison with SELEFA® Sense Blue Examination gloves, based on a life cycle analysis (LCA)
- ▶ High-performance, lightweight gloves that help reduce waste weight
- ▶ AQL 1.0, ensuring high quality
- ▶ Extended cuff for safer use
- ▶ Accelerator-free gloves with low levels of chemical residues, suitable for users with type IV allergy
- ▶ Latex-free and therefore suitable for those with latex / type I allergy
- ▶ Powder-free: reduces risk of allergy and minimises particulate contamination
- ▶ Suitable for contact with all type of foods
- ▶ Light, thin and strong



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
711005	Examination gloves, Nitrile PROCEED X-Long, XS, Green	XS	290 mm	100 / 1000
711006	Examination gloves, Nitrile PROCEED X-Long, S, Green	S	290 mm	100 / 1000
711007	Examination gloves, Nitrile PROCEED X-Long, M, Green	M	290 mm	100 / 1000
711008	Examination gloves, Nitrile PROCEED X-Long, L, Green	L	290 mm	100 / 1000
711009	Examination gloves, Nitrile PROCEED X-Long, XL, Green	XL	290 mm	100 / 1000

Sales Unit: Transport Carton



The Embra® Way

how we help build a positive future through responsible, futureproof healthcare products and solutions.



Embra products are designed to lower the environmental impact while heightening quality for patients and caregivers. This is The Embra® Way. It's our way of caring for patients, caregivers and the future. Ask your Onemed colleague for more information on how we together can strive to reduce our environmental impact and improve our sustainability performance.



Scan the QRcode and read more about Embra®





evercare[®] Examination Gloves, Nitrile SKY

REF: 6400, 6401, 6402, 6403, 6404

Nitrile examination glove and protective glove for single use, blue.

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ High functional resistance.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Excellent viral barrier properties, protects against a wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6400	Examination Gloves, Nitrile SKY, XS, Blue	XS	240 mm	100 / 1000
6401	Examination Gloves, Nitrile SKY, S, Blue	S	240 mm	100 / 1000
6402	Examination Gloves, Nitrile SKY, M, Blue	M	240 mm	100 / 1000
6403	Examination Gloves, Nitrile SKY, L, Blue	L	240 mm	100 / 1000
6404	Examination Gloves, Nitrile SKY, XL, Blue	XL	240 mm	100 / 1000

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile ROSA

REF: 6405, 6406, 6407, 6408, 6409

Nitrile examination glove and protective glove for single use, pink

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Very good viral barrier properties, protects against wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6405	Examination Gloves, Nitrile ROSA, XS, Pink	XS	240 mm	150 / 1500
6406	Examination Gloves, Nitrile ROSA, S, Pink	S	240 mm	150 / 1500
6407	Examination Gloves, Nitrile ROSA, M, Pink	M	240 mm	150 / 1500
6408	Examination Gloves, Nitrile ROSA, L, Pink	L	240 mm	150 / 1500
6409	Examination Gloves, Nitrile ROSA, XL, Pink	XL	240 mm	150 / 1500

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile SVEA

REF: 6415, 6416, 6417, 6418, 6419

Nitrile examination glove and protective glove for single use, blue/lilac.

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Well recognized examination glove color for ward.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Very good viral barrier properties, protects against a wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6415	Examination Gloves, Nitrile SVEA, XS, Blue/Lilac		240 mm	150 / 1500
6416	Examination Gloves, Nitrile SVEA, S, Blue/Lilac		240 mm	150 / 1500
6417	Examination Gloves, Nitrile SVEA, M, Blue/Lilac		240 mm	150 / 1500
6418	Examination Gloves, Nitrile SVEA, L, Blue/Lilac		240 mm	150 / 1500
6419	Examination Gloves, Nitrile SVEA, XL, Blue/Lilac		240 mm	150 / 1500

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile SAFE X-LONG

REF: 6420, 6421, 6422, 6423, 6424

Nitrile examination glove and protective glove for single use, x-long, blue.

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Strong, soft and flexible glove with excellent fit.
- ▶ Long cuff for extended protection of wrist.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Suitable for use in chemotherapy.
- ▶ Excellent viral barrier properties, protects against a wide variety of chemicals.
- ▶ The glove's durability has also been tested in cold temperature to confirm the suitability for use under cold working conditions outdoors such as e.g. ambulance work.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6420	Examination Gloves, Nitrile SAFE X-LONG, XS, Blue		295 mm	100 / 1000
6421	Examination Gloves, Nitrile SAFE X-LONG, S, Blue		295 mm	100 / 1000
6422	Examination Gloves, Nitrile SAFE X-LONG, M, Blue		295 mm	100 / 1000
6423	Examination Gloves, Nitrile SAFE X-LONG, L, Blue		295 mm	100 / 1000
6424	Examination Gloves, Nitrile SAFE X-LONG, XL, Blue		295 mm	100 / 1000

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile FROST

REF: 6425, 6426, 6427, 6428, 6429

Nitrile examination glove and protective glove for single use, white.

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Very good viral barrier properties, protects against wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6425	Examination Gloves, Nitrile FROST, XS, White		240 mm	150 / 1500
6426	Examination Gloves, Nitrile FROST, S, White		240 mm	150 / 1500
6427	Examination Gloves, Nitrile FROST, M, White		240 mm	150 / 1500
6428	Examination Gloves, Nitrile FROST, L, White		240 mm	150 / 1500
6429	Examination Gloves, Nitrile FROST, XL, White		240 mm	150 / 1500

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile LIME

REF: 6455, 6456, 6457, 6458, 6459


Nitrile examination glove and protective glove for single use, accelerator-free, green

Key Features

- ▶ Accelerator-free, powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Smooth surface at the fingertips for better feel.
- ▶ Highly elastic and smooth for excellent tactile sense
- ▶ Fits comfortably, easy to don and remove.
- ▶ Green color identifying of skin friendly material for use.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Good alternative to those suffering from type I allergies and accelerator-related type IV allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6455	Examination Gloves, Nitrile LIME, XS, Green		240 mm	150 / 1500
6456	Examination Gloves, Nitrile LIME, S, Green		240 mm	150 / 1500
6457	Examination Gloves, Nitrile LIME, M, Green		240 mm	150 / 1500
6458	Examination Gloves, Nitrile LIME, L, Green		240 mm	150 / 1500
6459	Examination Gloves, Nitrile LIME, XL, Green		240 mm	150 / 1500

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile SATIN

REF: 6490, 6491, 6492, 6493, 6494

Nitrile examination glove and protective glove for single use, white.

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ White color gives a signal of neutrality and purity.
- ▶ Very good viral barrier properties, protects against wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6490	Examination Gloves, Nitrile SATIN, XS, White		240 mm	150 / 1500
6491	Examination Gloves, Nitrile SATIN, S, White		240 mm	150 / 1500
6492	Examination Gloves, Nitrile SATIN, M, White		240 mm	150 / 1500
6493	Examination Gloves, Nitrile SATIN, L, White		240 mm	150 / 1500
6494	Examination Gloves, Nitrile SATIN, XL, White		240 mm	150 / 1500

Sales Unit: Transport Carton



evercare[®] Examination Gloves, Nitrile SATIN X-LONG

REF: 6520, 6521, 6522, 6523, 6524

Nitrile examination gloves and protective gloves for single use, x-long, white

Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Long cuff for extended protection of wrist.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Very good viral barrier properties, protects against wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



Recommendations for use



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
6520	Examination Gloves, Nitrile SATIN X-LONG, XS, White	XS	100 / 1000
6521	Examination Gloves, Nitrile SATIN X-LONG, S, White	S	100 / 1000
6522	Examination Gloves, Nitrile SATIN X-LONG, M, White	M	100 / 1000
6523	Examination Gloves, Nitrile SATIN X-LONG, L, White	L	100 / 1000
6524	Examination Gloves, Nitrile SATIN X-LONG, XL, White	XL	100 / 1000

Sales Unit: Transport Carton



SELEFA® Examination Gloves Nitrile SMART, White

REF: 210252XS, 210252S, 210252M, 210252L, 210252XL

Nitrile examination gloves and protective gloves for single use, white.

Key Features

- ▶ Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210252XS	Examination Gloves, Nitrile SMART, XS, White		240 mm	200 / 2000
210252S	Examination Gloves, Nitrile SMART, S, White		240 mm	200 / 2000
210252M	Examination Gloves, Nitrile SMART, M, White		240 mm	200 / 2000
210252L	Examination Gloves, Nitrile SMART, L, White		240 mm	200 / 2000
210252XL	Examination Gloves, Nitrile SMART, XL, White		240 mm	180 / 1800

Sales Unit: Transport Carton



SELEFA® Examination Gloves, Nitrile SMART, White

REF: 210256XS, 210256S, 210256M, 210256L, 210256XL

Nitrile examination gloves and protective gloves for single use, white.

Key Features

- ▶ Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210256XS	Examination Gloves, Nitrile SMART, XS, White		240 mm	150 / 1500
210256S	Examination Gloves, Nitrile SMART, S, White		240 mm	150 / 1500
210256M	Examination Gloves, Nitrile SMART, M, White		240 mm	150 / 1500
210256L	Examination Gloves, Nitrile SMART, L, White		240 mm	150 / 1500
210256XL	Examination Gloves, Nitrile SMART, XL, White		240 mm	150 / 1500

Sales Unit: Transport Carton



SELEFA® Examination Gloves, Nitrile SENSE, Blue

REF: 210258XS, 210258S, 210258M, 210258L, 210258XL

Nitrile examination gloves and protective gloves for single use, blue.

Key Features

- ▶ Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210258XS	Examination Gloves, Nitrile SENSE, XS, Blue		240 mm	200 / 2000
210258S	Examination Gloves, Nitrile SENSE, S, Blue		240 mm	200 / 2000
210258M	Examination Gloves, Nitrile SENSE, M, Blue		240 mm	200 / 2000
210258L	Examination Gloves, Nitrile SENSE, L, Blue		240 mm	200 / 2000
210258XL	Examination Gloves, Nitrile SENSE, XL, Blue		240 mm	180 / 1800

Sales Unit: Transport Carton



SELEFA® Examination Gloves Nitrile SENSE, Black

REF: 210259XS, 210259S, 210259M, 210259L, 210259XL


Nitrile examination gloves and protective gloves for single use, black.

Key Features

- ▶ Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210259XS	Examination Gloves, Nitrile SENSE, XS, Black		240 mm	200 / 2000
210259S	Examination Gloves, Nitrile SENSE, S, Black		240 mm	200 / 2000
210259M	Examination Gloves, Nitrile SENSE, M, Black		240 mm	200 / 2000
210259L	Examination Gloves, Nitrile SENSE, L, Black		240 mm	200 / 2000
210259XL	Examination Gloves, Nitrile SENSE, XL, Black		240 mm	180 / 1800

Sales Unit: Transport Carton



SELEFA Examination Gloves Nitrile SENSE LIGHT, White

REF: 210262XS, 210262S, 210262M, 210262L, 210262XL

Nitrile examination gloves and protective gloves for single use, white.

Key Features

- ▶ Thin nitrile glove allowing comfort and sensitivity to the user.
- ▶ Powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Good viral barrier properties, protects against chemicals.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Free of latex, powder and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210262XS	Examination gloves, Nitrile SENSE LIGHT, XS, White		240 mm	200 / 2000
210262S	Examination gloves, Nitrile SENSE LIGHT, S, White		240 mm	200 / 2000
210262M	Examination gloves, Nitrile SENSE LIGHT, M, White		240 mm	200 / 2000
210262L	Examination gloves, Nitrile SENSE LIGHT, L, White		240 mm	200 / 2000
210262XL	Examination gloves, Nitrile SENSE LIGHT, XL, White		240 mm	180 / 1800

Sales Unit: Transport Carton



SELEFA® Examination Gloves, Nitrile AQUA

REF: 6011055, 6011065, 6011075, 6011085, 6011095

Nitrile examination gloves and protective gloves for single use, blue

Key Features

- ▶ Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips.
- ▶ Fits comfortably, easy to don and remove.
- ▶ Good viral barrier properties, protects against variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.
- ▶ The gloves are suitable for all food contact.




Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	Pack Count dispenser/carton
6011055	Examination Gloves, Nitrile AQUA, XS, Blue		240 mm	100 / 1000
6011065	Examination Gloves, Nitrile AQUA, S, Blue		240 mm	100 / 1000
6011075	Examination Gloves, Nitrile AQUA, M, Blue		240 mm	100 / 1000
6011085	Examination Gloves, Nitrile AQUA, L, Blue		240 mm	100 / 1000
6011095	Examination Gloves, Nitrile AQUA, XL, Blue		240 mm	100 / 1000

Sales Unit: Transport Carton

Vinyl Exam Gloves



Product Data Sheet



evercare[®] Examination Gloves, Vinyl SENSITIVE
REF: 6350, 6351, 6352, 6353, 6354

Powder-free and phthalate-free vinyl examination glove for single use.

Key Features

- ▶ Ambidextrous with rolled rim for extra cuff strength
- ▶ Phthalate-free, powder-free
- ▶ Latex-free
- ▶ Good alternative to those suffering from type I or accelerator relating type IV allergy



Recommendations for use








Basic cleaning



Basic procedures
and wards



Wash

REF	DESCRIPTION	SIZE CODE	Overall Length	Pack Count dispenser/carton
6350	Examination Gloves, Vinyl SENSITIVE, XS, Clear		240 mm	100 / 1000
6351	Examination Gloves, Vinyl SENSITIVE, S, Clear		240 mm	100 / 1000
6352	Examination Gloves, Vinyl SENSITIVE, M, Clear		240 mm	100 / 1000
6353	Examination Gloves, Vinyl SENSITIVE, L, Clear		240 mm	100 / 1000
6354	Examination Gloves, Vinyl SENSITIVE, XL, Clear		240 mm	100 / 1000

Sales Unit: Transport Carton



SELEFA® Examination Gloves, Vinyl REGULAR, Phthalate-Free

REF: 210286XS, 210286S, 210286M, 210286L, 210286XL

Powder-free and phthalate-free vinyl examination glove for single use.

Key Features

- ▶ Ambidextrous and with rolled rim for extra cuff strength
- ▶ Phthalate-free, powder-free and latex-free
- ▶ Good alternative to those suffering from type I and accelerator related type IV allergies



Recommendations for use



REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
210286XS	Examination Gloves, Vinyl REGULAR, Phthalate-Free, XS, Clear		240 mm	100 / 1000
210286S	Examination Gloves, Vinyl REGULAR, Phthalate-Free, S, Clear		240 mm	100 / 1000
210286M	Examination Gloves, Vinyl REGULAR, Phthalate-Free, M, Clear		240 mm	100 / 1000
210286L	Examination Gloves, Vinyl REGULAR, Phthalate-Free, L, Clear		240 mm	100 / 1000
210286XL	Examination Gloves, Vinyl REGULAR, Phthalate-Free, XL, Clear		240 mm	100 / 1000

Sales Unit: Transport Carton



Embra® Latex Exam Glove



Embra® Examination Gloves, Latex FLEX, Powder-Free
REF: 821000, 821001, 821002, 821003, 821004

Disposable Examination Gloves, Latex, Powder-Free.

Key Features

- ▶ FSC-certified biobased material- Natural rubber Latex
- ▶ Elastic examination glove made of natural latex material with low latex protein content.
- ▶ Powder-free, ambidextrous with rolled rim and textured fingertips.
- ▶ Strong glove with a good fit and comfort.
- ▶ Excellent grip.
- ▶ Inner coated for easy donning.



Recommendations for use



Basic cleaning



Basic procedures
and wards

REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
821000	Examination Gloves, Latex FLEX, Powder-Free, XS, Off white		240 mm	100 / 1000
821001	Examination Gloves, Latex FLEX, Powder-Free, S, Off white		240 mm	100 / 1000
821002	Examination Gloves, Latex FLEX, Powder-Free, M, Off white		240 mm	100 / 1000
821003	Examination Gloves, Latex FLEX, Powder-Free, L, Off white		240 mm	100 / 1000
821004	Examination Gloves, Latex FLEX, Powder-Free, XL, Off white		240 mm	90 / 900

Sales Unit: Transport Carton

SELEFA® Copolymer Glove



SELEFA® Powder-Free Paper Backed Copolymer Gloves, Sterile REF: 6031865, 6031875, 6031885

A sterile, ambidextrous, powder & protein-free co-polymer glove attached to paper.

Key Features

- ▶ Paper backed sterile EVA Copolymer glove
- ▶ One glove per pouch



Recommendations for use




Basic procedures
and wards

REF	DESCRIPTION	SIZE	Overall Length	Pack Count dispenser/carton
6031865	Powder-Free Paper Backed Copolymer Glove, Sterile, Transparent		280 mm	100 / 2000
6031875	Powder-Free Paper Backed Copolymer Glove, Sterile, Transparent		280 mm	100 / 2000
6031885	Powder-Free Paper Backed Copolymer Glove, Sterile, Transparent		280 mm	100 / 2000

Sales Unit: Transport Carton

SELEFA® TPE Glove



SELEFA® Smooth Stretch TPE Gloves

REF: 210352S, 210352M, 210352L, 210352XL

Disposable Smooth Stretch TPE Gloves used primarily in food handling, general care, salons / spas, electronic assembly, laboratory work, packaging, automotive, janitorial, general pharmaceutical and home use.

Key Features

- ▶ Highly elastic, skin friendly, odourless and durable.
- ▶ Thermoplastic elastomer (TPE).
- ▶ Ambidextrous, with good fit on hand and fingertips.
- ▶ Easy to don and take off.
- ▶ Smooth surface and non-beaded cuff.
- ▶ Suitable for all type of food handling.



Recommendations for use



REF	DESCRIPTION	SIZE CODE	Overall Length	PACK COUNT dispenser/carton
210352S	Smooth Stretch TPE Gloves, S, Transparent		260 mm	200 / 2000
210352M	Smooth Stretch TPE Gloves, M, Transparent		260 mm	200 / 2000
210352L	Smooth Stretch TPE Gloves, L, Transparent		260 mm	200 / 2000
210352XL	Smooth Stretch TPE Gloves, XL, Transparent		260 mm	200 / 2000

Sales Unit: Transport Carton

evercare[®] medical


We care more



www.evercare-medical.com




Follow us on **LinkedIn**



Read latest **news**

Embra [®]
Caring for the future

 **evercare** [®]
Solutions for premium care

 **SELEFA** [®]
Affordable solutions for care

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.